

Real Time PCR Mustard DNA Detection Kit

Test system for the qualitative detection of mustard DNA in food products by PCR Real time

Product code: IC-02-1036 (50 tests) / IC-02-1038 (25 tests)

Brief description

MustardKit Real Time PCR provides reagents for the qualitative detection of mustard DNA in several food products, fresh and processed. The PCR Real Time kit amplifies a DNA fragment that is present solely in mustard. The amplified DNA segment is detected by hybridisation with a probe labelled with fluorescent dyes. The increase in fluorescence is continuously measured in a PCR real-time detection instrument.

Mustard is considered allergenic food and it is explicitly mentioned in the European Food Labelling Directive.

MustardKit Real Time PCR is an useful tool to monitor the food allergens to ensure compliance with food labelling (2003/89/EC Directive) and to improve consumer protection.




Figure 1: Test sensitivity for Mustardkit DNA Detection kit.

— 1000 copies mustard genome
 — 100 copies mustard genome
 — 10 copies mustard genome
 — 1 copy mustard genome

Technical features

Number of tests	25/50 target DNA specific reactions
Kit components	Mix <i>Test Mustard</i> (with duplex inhibition control) – DNA positive control – sterile H ₂ O DNase free
Specificity	the kit has been tested with DNA extracted from both raw material and processed food. The specificity of the system has been validated with several other species normally used in food production (see the Panel below). No cross-reactivity effects have been revealed.
Limit of Detection	1 copy of mustard haploid genome, equal to roughly 0.5 pg of <i>Sinapis alba</i> DNA (see Fig.1).
Detection	Probe labelled with fluorescent dyes Taqman® - FAM and JOE

Specificity control group tested **negative** for cross-reactivity:

Almond (<i>Prunus dulcis</i>)	Celery (<i>Apium graveolens</i>)	Peanut (<i>Arachis hypogaea</i>)
Barley (<i>Hordeum vulgare</i>)	Chicken (<i>Gallus gallus</i>)	Pistachio nut (<i>Pistacia vera</i>)
Black tiger prawn (<i>Penaeus monodon</i>)	Cod (<i>Merluccius capensis</i>)	Radish (<i>Raphanus sativus</i>)
Bovine (<i>Bos taurus</i>)	Common wheat (<i>Triticum aestivus</i>)	Rucola (<i>Eruca sativa</i>)
Brasilian walnut (<i>Bertholletia excelsa</i>)	<i>aestivus</i>)	Rye (<i>Secale cereale</i>)
Broccoli (<i>Brassica oleracea convar.</i>)	Hazelnut (<i>Coryllus avellana</i>)	Savoy cabbage (<i>Brassica oleracea</i>)
<i>Botrytis cult. Italica</i>)	Lupin (<i>Lupinus luteus</i>)	<i>capitata conv cult. Savoy</i>)
Cauliflower (<i>Brassica oleracea convar.</i>)	Man (<i>Homo sapiens</i>)	Sesame (<i>Sesamum indicum</i>)
<i>Botrytis</i>)	Mussel (<i>Mytilus chilensis</i>)	Soybean (<i>Glycine max</i>)
Cashew nut (<i>Anacardium occidentale</i>)	Papaya (<i>Carica papaya</i>)	Walnut (<i>Juglans regia</i>)