

# WHEAT<sub>a</sub> KIT

## REAL TIME PCR DNA DETECTION KIT

### Test system for the qualitative and quantitative detection of *Triticum aestivum* in *Triticum turgidum* DNA by means of PCR Real time


Product code: IC-02-1046 (100 reactions)

#### Brief description

Cereal composition is a key factor in the quality and safety of food and feed. In particular, *durum* wheat is considered the election wheat to make pasta, given its unique properties like relatively high yellow pigment content and high protein content favourable for good cooking quality. The dough made from *durum* wheat has rheological properties ideally suited to the pasta manufacturing process.

Because of the possibility of accidental contamination occurring during either wheat harvest or storage and transport of grains and semolina, pasta is only officially regarded as impure when the common wheat level exceeds 3%, as stipulated by the European Commission regulation (1222/94, EC 1994).

**WHEAT<sub>a</sub>KIT DNA Detection Kit** allows the quantitative detection in samples like flour, semolina and pasta of common wheat "*Triticum aestivum*" DNA in *durum* wheat "*Triticum turgidum*" DNA.



#### Technical features

<b>Number of tests</b>	100 target DNA specific reactions.
<b>Kit content</b>	Mix <i>Triticum</i> – Standard S1, S2, S3, S4, S5 - sterile H <sub>2</sub> O DNase free.
<b>Specificity</b>	The test system is based on selective identification of a common genetic marker of <i>Triticum</i> genus (see Table 2 and Figure 2) and of a unique genetic marker of <i>Triticum aestivum</i> cultivars (see Table 1 and Figure 1).
<b>Sensibility</b>	up to 0.2% of <i>Triticum aestivum</i> in <i>Triticum turgidum</i> .
<b>Limit of Detection</b>	1 copy of haploid genome of <i>Triticum aestivum</i> .
<b>Detection</b>	probe labelled with fluorescent dyes Taqman® - FAM and JOE

#### Tab. 1 Species verified for the absence of amplification for the unique genetic marker of *Triticum aestivum* cultivars:

Barley (*Hordeum vulgare*), Buckwheat (*Fagopyrum esculentum*), Chickpea (*Cicer arietinum*), Corn (*Zea mays*), Durum wheat (*Triticum turgidum durum*), Festuca arundinacea (*Festuca urundinacea*), Flax (*Linum usatissimum*), Kamut (*Triticum turgidum polonicum*), Lentil (*Lens esculenta*), Millet (*Sorghum bicolor*), Oat (*Avena sativa*), Orchard grass (*Dactylis glomerata*), Perennial Ryegrass (*Lolium perenne*), Rapeseed (*Brassica napus*), Reddish purple (*Amaranthus spp*), Red mile (*Panicum miliaceum*), Rice (*Oryza sativa*), Rye (*Secale cereale*), Spelt (*Triticum turgidum dicoccum*), Sunflower (*Helianthus annuus*), Sugar cane (*Saccharum officinarum*), Soy (*Glycine max*), Triticale (X *Triticosecale*).

#### Tab. 2 Species verified for the amplification of common genetic marker of *Triticum* genus:

Common wheat (*Triticum aestivum*), Durum wheat (*Triticum turgidum durum*), Kamut (*Triticum turgidum polonicum*), Spelt (*Triticum turgidum dicoccum*), Triticale (X *Triticosecale*).