

Real Time PCR Hazelnut DNA Detection Kit

Test system for the qualitative detection of hazelnut DNA in food products by PCR Real time

Product code: IC-02-1031 (50 tests) / IC-02-1033 (25 tests)

Brief description

HazelnutKit Real Time PCR provides reagents for the qualitative detection of hazelnut DNA in several food products, fresh and processed. The PCR Real Time kit amplifies a DNA fragment that is present solely in hazelnut. The amplified DNA segment is detected by hybridisation with a probe labelled with fluorescent dyes. The increase in fluorescence is continuously measured in a PCR real-time detection instrument.

Hazelnut is considered allergenic food and it is explicitly mentioned in the European Food Labelling Directive.

HazelnutKit Real Time PCR is an useful tool to monitor the food allergens to ensure compliance with food labelling (2003/89/EC Directive) and to improve consumer protection.





Figure 1: Test sensitivity for *Hazelnutkit* DNA Detection kit.

— 1000 genome copies *Corylus avellana*
— 100 genome copies *Corylus avellana*
— 10 genome copies *Corylus avellana*
— 1 genome copy *Corylus avellana*

Technical features

Number of tests	25/50 target DNA specific reactions
Kit components	Mix <i>Test Hazelnut</i> – Mix Test Inhibition – DNA positive control – sterile H ₂ O DNase free
Specificity	the kit has been tested with DNA extracted from both raw material and processed food. The specificity of the system has been validated with several other species normally used in food production (see the Panel below). No cross-reactivity effects have been revealed.
Limit of Detection	1 copy of hazelnut haploid genome, equal to roughly 0.48 pg of <i>Corylus avellana</i> DNA (see Fig.1).
Detection	Probe labelled with fluorescent dyes Taqman® - FAM

Specificity control group tested **negative** for cross-reactivity:

Almond (<i>Prunus dulcis</i>)	Currant (<i>Rubus vulgare</i>)	Pecan nut (<i>Carya illinoensis</i>)
Apple (<i>Malus domestica</i>)	Giuggiola (<i>Zizyphus jujuba</i>)	Peppermint (<i>Mentha piperita</i>)
Apricot (<i>Prunus armeniaca</i>)	Green Pea (<i>Pisum sativum</i>)	Pine nut (<i>Pinus pinea</i>)
Barley (<i>Hordeum vulgare</i>)	Lemon (<i>Citrus limon</i>)	Pistachio (<i>Pistacia vera</i>)
Bean (<i>Phaseolus vulgaris</i>)	Lentils (<i>Lens esculenta</i>)	Plum (<i>Prunus domestica</i>)
Blackberry (<i>Rubus fruticosus</i>)	Macadamia nut (<i>Macadamia integrifolia</i>)	Potato (<i>Solanum tuberosum</i>)
Black cherry (<i>Prunus cerasus</i>)	Maize (<i>Zea mays</i>)	Prickly pear (<i>Opuntia ficus indica</i>)
Bovine (<i>Bos taurus</i>)	Man (<i>Homo sapiens</i>)	Pumpkin (<i>Cucurbita pepo</i>)
Brazil nut (<i>Bertholletia excelsa</i>)	Mussel (<i>Mytilus chilensis</i>)	Raspberry (<i>Rubus idaeus</i>)
Broad bean (<i>Vicia faba</i>)	Mile (<i>Panicum miliaceum</i>)	Rice (<i>Oryza sativa</i>)
Buckwheat (<i>Fagopyrum esculentum</i>)	Mushroom (<i>Agaricus bisporus</i>)	Rye (<i>Secale cereale</i>)
Carrot (<i>Daucus carota</i>)	Mustard (<i>Brassica alba</i>)	Sesame (<i>Sesamum indicum</i>)
Cashew nut (<i>Anacardium occidentale</i>)	Nutmeg (<i>Myristica fragrans</i>)	Shrimp (<i>Peneaus monodon</i>)
Cherry (<i>Prunus avium</i>)	Oat (<i>Avena sativa</i>)	Soya (<i>Glycine max</i>)
Chestnut (<i>Castanea sativa</i>)	Orange (<i>Citrus sinensis</i>)	Spelt (<i>Triticum dicoccum</i>)
Chicken (<i>Gallus gallus</i>)	Parsley (<i>Petroselinum crispum</i>)	Triticale (<i>Triticosecalle</i>)
Celery (<i>Apium graveolens</i>)	Peach (<i>Prunus persica</i>)	Walnut (<i>Juglans regia</i>)
Coconut (<i>Cocos nucifera</i>)	Peanut (<i>Arachis hypogaea</i>)	Wheat (<i>Triticum aestivum</i>)
Cod (<i>Merluccius capensis</i>)	Pear (<i>Pyrus communis</i>)	White Lupin (<i>Lupinus albus</i>)
		Yeast (<i>Saccharomyces cerevisiae</i>)