

Buffer B Milk 10X

Extraction Buffer for samples containing poly-phenols

Code: IC-02-1112 / 15 samples

IC-02-1123 / 50 samples

Milk II Alert ELISA Kit



Description

Buffer B Milk allows the extraction of proteins from poly-phenols containing food matrices.

Required for the sample preparation in **Milk II Alert ELISA kit**.

Features

Content	Buffer B Milk 10X - 1 x 15 ml / 1 x 50 ml
Number of samples	15
Samples	<u>red fruits</u> : blueberries, currants, raspberries, strawberries, blackberries <u>spices</u> : mustard, chilly pepper, basil, curry <u>chocolate</u> <u>cocoa</u> and cocoa containing foodstuffs <u>nuts</u>
Storage	2 – 8°C
Stability	1 year
Related products	Milk II Alert ELISA kit for the detection of milk proteins in foodstuff (IC-02-1102 / IC-02-1103)
Instructions	the sample preparation is described in the user manual of Milk II Alert ELISA kit .