

EPA Extraction Powder

Extraction Powder Additive for proteins extraction from foodstuff containing:

- ▶ cocoa, in association with **SoyAlert ELISA kit**
- ▶ cocoa and red fruits in association with **HazelnutAlert ELISA kit**

Code: IC-02-1122

Description

EPA Extraction Powder allows the extraction of proteins from food matrices with high content of **polyphenols**, such cocoa and red fruits.

Required for the sample preparation in **SoyAlert ELISA kit** and **HazelnutAlert ELISA kit**

Features

Content	1 x 15 g
Number of samples	50 samples with SoyAlert ELISA kit 30 samples with HazelnutAlert ELISA kit
Samples	foodstuff containing polyphenols
Storage	2-8 °C
Stability	3 years
Related products	SoyAlert ELISA kit for the detection of soy proteins in foodstuff containing cocoa (IC-02-1125 / IC-02-1126) HazelnutAlert ELISA kit for the detection of hazelnut proteins in foodstuff containing cocoa and red fruits (IC-02-1117/ IC-02-1118)
Instructions	the sample preparation is described in related user manuals