

# HM Extraction Buffer

## Milk II Alert ELISA Kit

Extraction Buffer for protein extraction from matrices like **ham**

Code: IC-02-1109



### Description

**HM Extraction Buffer** allows the extraction of proteins from food matrices like **ham**.

Required for the sample preparation in **Milk II Alert ELISA kit**.

### Features

<b>Content</b>	<b>HM Buffer</b> ready to use - 1 x 50 ml.
<b>Number of samples</b>	50
<b>Samples</b>	Food matrices like ham.
<b>Storage</b>	Room Temperature.
<b>Stability</b>	2 years.
<b>Related products</b>	<b>Milk II Alert ELISA kit</b> for the detection of milk proteins in foodstuff (IC-02-1102 / IC-02-1103).
<b>Instructions</b>	the sample preparation is described in the user manual of <b>Milk II Alert ELISA kit</b> .