

WM Extraction powder

Milk II Alert ELISA Kit

Extraction Powder for protein extraction from matrices like **wines** and **vinegars**

Code: IC-02-1104



Description

WM Extraction Powder allows the extraction of proteins from food matrices like wines and vinegars.

Required for the sample preparation in **Milk II Alert ELISA kit**.

Features

Content	WM Powder - 1 x 25 g
Number of samples	25
Samples	<u>wines</u> : white and red <u>vinegars</u> : from white and red wines, apples
Storage	2 – 8°C
Stability	3 years
Related products	Milk II Alert ELISA kit for the detection of milk proteins in foodstuff (IC-02-1102 / IC-02-1103)
Instructions	the sample preparation is described in the user manual of Milk II Alert ELISA kit .